



The Dolomite Centre Ltd
 Unit 1, Anglian Business Park, Orchard Road,
 Royston, Hertfordshire, SG8 5TW, UK

T: +44 (0)1763 242491
F: +44 (0)1763 246125
E: sales@dolomite-microfluidics.com
W: www.dolomite-microfluidics.com

Dolomite Microfluidics – North America Office
 Blacktrace Inc, 29 Albion Place
 Charlestown, MA 02129, USA


C: 617 803 6655
T: 617 848 1121
F: 617 500 0136
E: salesus@dolomite-microfluidics.com

INFORMATION SHEET


| | | | |
|-----------|-----------------------|-------------|---------|
| Part name | Mitos micromixer chip | Part number | 3000144 |
|-----------|-----------------------|-------------|---------|

Description

The Mitos micromixer chip is a glass microfluidic device designed for rapid mixing of two or three fluid streams. Fast mixing times are required for the study of reaction kinetics, sample dilution, improving reaction selectivity, rapid crystallisation and nanoparticle synthesis. There are two independent micromixer channels on the chip and two headers allowing connection with the Mitos 4-way edge connector (3000024).



Left: Mitos micromixer chip (3000144)
 Below: Micromixer chip with Edge Connectors



Benefits


- Extremely rapid mixing across a range of flow rates
- Low dead volume
- High visibility (excellent access for optics)
- Quick connect / disconnect
- Wide temperature and pressure range
- Excellent chemical compatibility

| | Specification | Value |
|----|--|--|
| 1 | Number of inputs | 3 |
| 2 | Number of outputs | 1 |
| 3 | Internal channel cross section | 125 x 350 μm and 50 x 125 μm (depth x width) |
| 4 | Internal volume of micromixer | 8 μl |
| 5 | Internal volume of a mixing stage | 0.35 μl |
| 6 | Back pressure with 100 $\mu\text{l}/\text{min}$ flow (water) | 0.1 bar |
| 7 | Outside diameter of connection tubing | 1.6 mm (1/16 inch) |
| 8 | Inside diameter of connection tubing | 0.25 mm provided as standard |
| 9 | Connection tubing material | PTFE, FEP |
| 10 | Surface roughness of channels (R_a) | 5 nm |
| 11 | Chip size | 45 mm x 15 mm |
| 12 | Chip top layer thickness | 2 mm |
| 13 | Chip base layer thickness | 2 mm |
| 14 | Operating pressure | 30 Bar with edge connector |

| | | |
|----|-----------------------|---|
| 15 | Operating temperature | 250 °C with edge connector (the glass can be locally heated to higher temperatures) |
| 16 | Material | Glass |
| 17 | Fabrication process | HF etching and thermal bonding |

Micromixer performance

To measure mixing performance a transparent phenolphthalein pH marker was mixed with an alkali solution to produce a bright pink colour. This is shown in the image below.




The mixing time was measured at various total flow rates, as shown in the table below. The volumetric flow rate ratio between the two input streams was 1:1.

| Flow rate (microliters/min) | 5 | 10 | 20 | 40 | 80 | 160 | 320 | 640 | 1200 | 2000 | 3000 | 4000 | 5000 |
|-----------------------------|------|------|------|------|------|-----|-----|-----|------|------|------|------|------|
| Number of stages for mixing | 1 | 1 | 2 | 4 | 6 | 7 | 7 | 4 | 4 | 3 | 3 | 2 | 2 |
| Mixing time (milliseconds) | 4200 | 2100 | 2100 | 2100 | 1575 | 919 | 459 | 131 | 70 | 32 | 21 | 11 | 8 |



Micromixer chip channel layout

The channel geometry of the Micromixer chip is shown below:



Mixing mechanism

The Mitos micromixer is a static mixer (no moving parts), which at low flow rates creates lamination of the flow streams as shown in the diagram below:



The lamination of the streams reduces diffusion distances and hence improves mixing time. At high flow rates swirling occurs in the flow streams, reducing mixing time further.